

SL NO	TEST CODE	INVESTIGATION NAME	NOTE	TEMPERATURE	REPORTING TIME	COST
IMMUNOCHEMISTRY/HORMONES						
1.	T3	T3: TRIIODOTHYRONINE TOTAL	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
2.	T4	T4: TETRAIODOTHYRONINE TOTAL	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
3.	FT3-1	FT3: FREE TRIIODOTHYRONINE	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
4.	FT4-1	FT4: FREE TETRAIODOTHYRONINE	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
5.	TSH	TSH: THYROID STIMULATING HORMONE	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
6.	ATPO	ATPO: ANTI THYROPEROXIDASE ANTIBODY	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
7.	TFT	THYROID PROFILE: THYROID FUNCTION TEST (TFT- BASIC)	COLLECT 3.5 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1.5 ML SERUM.	2-8 ⁰ C	SAME DAY	
		T3, T4, TSH				
8.	TFTEXT	FT3: FREE TRIIODOTHYRONINE, FT4: FREE TETRAIODOTHYRONINE, ATPO & TFT BASIC	COLLECT 3.5 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1.5 ML SERUM.	2-8 ⁰ C	SAME DAY	
9.	FSH/LH/PRL	FSH: FOLLICLE-STIMULATING HORMONE, LH: LUTEINIZING HORMONE & PRL: PROLACTIN	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM. SPECIFY AGE, SEX AND DAY OF MENSTRUAL CYCLE ON TEST REQUISITION FORM	2-8 ⁰ C	SAME DAY	
10.	FSH	FSH: FOLLICLE-STIMULATING HORMONE	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM. SPECIFY AGE, SEX AND DAY OF MENSTRUAL CYCLE ON TEST REQUISITION FORM	2-8 ⁰ C	SAME DAY	
11.	LH	LH: LUTEINIZING HORMONE	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
12.	PRL	PRL: PROLACTIN	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	

13.	TEST	TESTOSTERON TOTAL	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
14.	CORTI	CORTISOL BASELINE	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
15.	CORTAM	CORTISOL MORNING	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
16.	CORTPM	CORTISOL EVENING	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
17.	CORTISOL1	CORTISOL 1 HR AFTER ACTH	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
18.	PROG	PROGESTERON	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
19.	INSULIN	INSULIN FASTING	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM. OVERNIGHT FASTING IS MANDATORY.	2-8 ⁰ C	SAME DAY	
20.	INSULIN1	INSULIN PP	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM. SAMPLES TO BE COLLECTED 2 HOURS POST MEAL.	2-8 ⁰ C	SAME DAY	
21.	INSULINR	INSULIN RANDOM	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
22.	BHCG	BETA HCG	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
23.	FNBHCG	FREE BETA HCG	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
24.	PSA	PSA TOTAL	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
25.	PAPP-A	PAPP-A	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
26.	CPEP	C-PEPTIDE	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
27.	IGE	IGE: IMMUNOGLOBULIN E TOTAL	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
28.	FER	FERRITIN	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
29.	AFP	AFP: ALPHA FETOPROTEIN	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
30.	TM	TRIPPLE MARKER	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	4 DAYS	
31.	DM	DOUBLE MARKER	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	4 DAYS	

32.	CA125	CA 125: OVARIAN CANCER MARKER	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 ⁰ C	SAME DAY	
-----	-------	-------------------------------	--	--------------------	----------	--